

PF Series – Multipurpose Examination lights

This high intensity general exam light is flexible, durable and easy to use. Ideal for general purpose examinations.

PF FEATURES

- >>> Swivel joint lamp head allows for easy vertical and horizontal adjustment of light
- >>> Adjustable arm levers allow for easy positioning
- >>> Choice of articulated arm or flexible gooseneck models
- >>> Available in ceiling, wall, direct mount, and mobile configurations

TECHNICAL DATA

Power supply: 120V 60 Hz AC or 230-240V 50 Hz AC

Light intensity: Up to 4,000 Footcandles/43,200 Lux

Light source: Halogen bulb – Dichroic 12V/50W/10°

Light color temp: Up to 4,000 Kelvin

Color: White – RAL 9002

Configurations:	Wall	Ceiling	Direct Mount	Mobile
Weight:	17 lbs	17 lbs	4.5 lbs	23 lbs

Wall	Part No.: PF2015-44
Wall (not shown)	Part No.: PF2015-44PI
Ceiling	Part No.: PF2017-44
Direct Mount (not shown)	Part No.: PF28010
Mobile , gooseneck arm	Part No.: PF2006-32

