

## PF Series – Multipurpose Examination lights

This high intensity general exam light is flexible, durable and easy to use. Ideal for general purpose examinations.

### PF FEATURES

- >>> Swivel joint lamp head allows for easy vertical and horizontal adjustment of light
- >>> Adjustable arm levers allow for easy positioning
- >>> Choice of articulated arm or flexible gooseneck models
- >>> Available in ceiling, wall, direct mount, and mobile configurations

### TECHNICAL DATA

**Power supply:** 120V 60 Hz AC or 230-240V 50 Hz AC

**Light intensity:** Up to 4,000 Footcandles/43,200 Lux

**Light source:** Halogen bulb – Dichroic 12V/50W/10°

**Light color temp:** Up to 4,000 Kelvin

**Color:** White – RAL 9002

Configurations:	Wall	Ceiling	Direct Mount	Mobile
<b>Weight:</b>	17 lbs	17 lbs	4.5 lbs	23 lbs

<b>Wall</b>	Part No.: <b>PF2015-44</b>
<b>Wall</b> (not shown)	Part No.: <b>PF2015-44PI</b>
<b>Ceiling</b>	Part No.: <b>PF2017-44</b>
<b>Direct Mount</b> (not shown)	Part No.: <b>PF28010</b>
<b>Mobile</b> , gooseneck arm	Part No.: <b>PF2006-32</b>

