

## SoloSpot™

The SoloSpot is ideal for applications where a focused, spot light is required. The specially designed lamp head offers cool light and a focused, intense spot. This light is perfect for direct attachment to medical equipment such as patient tables and chairs. Ideal for OB-GYN, dental, and surgical applications.

### SOLOSPOT™ FEATURES

- >>> Unique, double filtered lens and color correction filter with lampshade extension help keep light cool on skin and illuminate with natural skin tones to improve tissue clarity
- >>> Excellent maneuverability allows for easy adjustment of light
- >>> Intense 20W Halogen light source (average 2,000 hours)
- >>> 3.5-inch diameter spot pattern @ 20-inches
- >>> Available in wall, ceiling, c-clamp and mobile configurations



**Mobile**

### TECHNICAL DATA

**Power supply:** 120V 60 Hz AC or 230-240V 50 Hz AC  
**Light intensity:** Up to 3,000 Footcandles/32,400 Lux  
**Light source:** Halogen bulb – 12V/20W  
**Light color temp:** 4,000 Kelvin  
**Color:** White – RAL 9002

Configurations:	Wall	Ceiling	C-clamp	Mobile
<b>Weight:</b>	6 lbs	16 lbs	6 lbs	15 lbs

<b>Wall</b>	Part No.: <b>HM2005-37EXT</b>
<b>C-clamp</b>	Part No.: <b>HM2002-38EXT</b>
<b>Mobile</b>	Part No.: <b>HM2006-38EXT</b>
<b>Ceiling</b>	Part No.: <b>HM2017-44EXT</b>



**C-Clamp**



**Wall**