


PF Series MRI Exam

This unique exam light is designed for use in any MRI facility. Made of durable, MRI compatible materials, this light is ideal for imaging centers, radiology departments, and CAT/PET Scan labs and clinics.

MRI EXAM FEATURES

- >>> Tested for MRI compatibility
- >>> Mobile caster base with easy-glide wheels
- >>> Powerful, 50W Halogen light source
- >>> Adjustable jointed arm for easy positioning
- >>> Weighted base for added stability



TECHNICAL DATA

Power supply: 120V 60 Hz AC or 230-240V 50 Hz AC
Light intensity: up to 4,300 Footcandles/46,000 Lux
Light source: Halogen bulb – Dichroic 12V/50W/10°
Light color temp: Up to 4,000 Kelvin
Color: White – RAL 9002

Configurations:	Mobile
Weight:	40 lbs

Mobile	Part No.: PF28016-M
---------------	----------------------------