


TriStar™

Excellent performance and great value is what the TriStar™ has to offer. Ideal for medical procedures in environments where space is limited, and multiple procedures are performed. The unique telescoping pole offers height adjustment in the mobile configuration and is one of the only products in the market with this functionality. The wall model offers extensive reach with a unique extension arm design.

TRISTAR™ FEATURES

- >>> Unique mobile stand, raises from 27 to 80 inches via one hand control for quick height adjustment
- >>> Casters move easily across all surfaces; base design prevents tipping
- >>> Wall model offers 30 inch extension arm, providing 55 inches of reach to beam center
- >>> (3) 20W long life Halogen bulbs (average over 4,000 hours)
- >>> Available in wall and mobile configurations



TECHNICAL DATA

Power supply: 120V 60 Hz AC or 230-240V 50 Hz AC

Light intensity: Up to 4,000 Footcandles/43,200 Lux

Light source: Halogen bulb – Dichroic 12V/3x20W/10°

Light color temperature: Up to 4,100 Kelvin

Medical standards: FDA Class II

Color: White – RAL 9002

Configurations:	Wall	Mobile
------------------------	------	--------

Weight:	8 lbs	25 lbs
----------------	-------	--------

Wall	Part No.: TS2015W
-------------	--------------------------

Mobile	Part No.: TS2015M
---------------	--------------------------

Extra Sterilizable Handle	Part No.: 18544
----------------------------------	------------------------

Color Corrective Filter	Part No.: 18676
--------------------------------	------------------------

